

## Comfort Series Material List All Parts are LATEX FREE

Product	Material
<b>ComfortGel</b>	
Outer Flap	Silicone / Polycarbonate
Gel Cushion Outer Shell	Polyurethane Film
Cushion Gel	Polyurethane Gel / EVA
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Elbow Split Washer	Polycarbonate and/or PETG
Forehead Support Arm and Bracket	Polycarbonate
Forehead Cushion	Silicone
Swivel Headgear Clip	Acetal
Deluxe Headgear	Velstretch / Lycra Laminate
<b>ComfortSelect</b>	
Cushion	Silicone Rubber
White Port Cap	Flexible Vinyl (PVC)
Split Washer	Polycarbonate and/or PETG
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Forehead Cushion	Silicone
Deluxe Headgear	Velstretch / Lycra Laminate
<b>ComfortLite</b>	
Cushion (Simple or DirectSeal)	Silicone
Simple cushion elbow	Polycarbonate
DirectSeal Cushion Cradle	Polycarbonate
DirectSeal Cushion Support Band	Aluminum
Crown Swivel	Polycarbonate
Flex Tube	EVA copolymer
Pivot Tube Insert	Polycarbonate
Height Adjustment Nut	Acetal
Height Adjustment Tube	Polycarbonate
Pivot Mount Strap	Polycarbonate
Pivot Tube	Polypropylene
Headband and Crown Strap	High Density Polyethylene
Headband and Crown Strap Liner	Micromesh Tricot with Akwadyne, Micro Ban foam
Rear Strap	Velcro Hook and Loop
Rear Strap Clip	High Density Polyethylene
<b>ComfortFull</b>	
Cushion	Silicone
Cushion Ring	PPE/PA Blend
Faceplate	Polycarbonate
Valve Body	Polycarbonate
Forehead Cushion	Silicone
ComfortHeadgear	Nylon / Lycra Laminate

## Additional Repironics Masks & Accessories Material List All Parts are LATEX FREE

Product	Material
<b>ComfortClassic</b>	
Cushion	Silicone
Faceplate	Polycarbonate
Split Washer	Polycarbonate and/or PETG
White Port Cap	Flexible Vinyl (PVC)
Swivel Elbow w/Exhalation Port	Polycarbonate
GEL Forehead Spacer	Polyurethane Elastic Film & GEL
Deluxe Headgear	Velstretch / Lycra Laminate
<b>Profile Lite / GoldSeal</b>	
Cushion Outer Shell	Polyurethane Elastic Film
Cushion Gel	Polyurethane Gel
White Port Cap	Flexible Vinyl (PVC)
Split Elbow	Polycarbonate and/or PETG
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Spacer	White Polyethyl Foam
Deluxe Headgear	Velstretch / Lycra Laminate
<b>Simplicity Nasal Mask</b>	
Faceplate	Polycarbonate
Tubing	Ethylene Vinyl Acetate (EVA)
Cushion	Silicone
Simplicity Headgear	Neoprene/Nylon
Deluxe Softcap	Velstretch / Lycra-Nylon
<b>Total Face Mask</b>	
Cushion	Silicone
Faceplate	Polycarbonate
Seal Ring	Alcryn
Valve Body	Polycarbonate
Total Headgear	Neoprene/Nylon
<b>Spectrum Reusable Full Face</b>	
Cushion	Polyurethane
Faceplate	Propionate
Seal Ring	Kraton Elastomer
Valve Body	Polycarbonate
Forehead Pad	White Polyethyl Foam
<b>Softcaps</b>	
Cap	Lycra
Strap Loop / Hook	Velstretch / Nylon Hook
<b>Simple Strap</b>	
Mesh	Polyester
Strap Loop / Hook	Velstretch / Nylon Hook